



# Recombinant Citrate lyase subunit beta (citE)

<b>Product Code</b>	CSB-EP364335EGX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A9I2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MISASLQQRK TRTRRSMLFV PGANAAMVSN SFIYPADALM FDLEDSVALR EKDTARRMVY HALQHPLYRD IETIVRVNAL DSEWGVNDLE AVVRGGADV RLPKTDTAQD VLDIEKEILR IEKACGREPG STGLLAAIES PLGITRAVEI AHASERLIGI ALGAEDYVRN LRTERSPEGT ELLFARCSIL QAARSAGIQA FDTVYSDANN EAGFLQEAAH IKQLGFDGKS LINPRQIDLL HNLYAPTQKE VDHARRVVEA AEAAAREGLG VVSLNGKMVD GPVIDRARLV LSRAELSGIR EE
<b>Source</b>	E.coli
<b>Target Names</b>	citE
<b>Protein Names</b>	Recommended name: Citrate lyase subunit beta Short name= Citrase beta chain EC= 4.1.3.6 Alternative name(s): Citrate (pro-3S)-lyase subunit beta Citryl-CoA lyase subunit EC= 4.1.3.34
<b>Expression Region</b>	1-302
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.