



# Recombinant Escherichia coli Ribonuclease G (rng)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP364339ENV  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | P0A9J0   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Escherichia coli (strain K12)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | TAELLVNVT PSETRVAYID GGILQEIHIE REARRGIVGN IYKGRVSRVL<br>PGMQAAFVDI GLDKAAFLHA SDIMPHTECV AGEEQKQFTV RDISELVRQG<br>QDLMVQVVKD PLGTKGARLT TDITLPSRYL VFMPGASHVG VSQRIESESE<br>RERLKKVVAE YCDEQGGFII RTAAEGVGEA ELASDAAYLK RVWTKVMERK<br>KRPQTRYQLY GELALAQRVL RDFADAELDR IRVDSRLTYE ALLEFTSEYI<br>PEMTSKLEHY TGRQPIFDLF DVENEIQRAL ERKVELKSGG YLIIDQTEAM<br>TTVDINTGAF VGHRNLDDTI FNTNIEATQA IARQLRLRNL GGIIDFID<br>MNNEDHRRRV LHSLEQALS DRVKTSVNGF SALGLVEMTR KRTRESIEHV<br>LCNECPTCHG RGTVKTIVETV CYEIMREIVR VHHAYDSDRF LUYASPAVAE<br>ALKGEESHSL AEVEIFVGKQ VKVQIEPLYN QEQFDVMM |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | rng  |
| <b>Protein Names</b>     | Recommended name: Ribonuclease G Short name= RNase G EC= 3.1.26.-<br>Alternative name(s): Cytoplasmic axial filament protein   |
| <b>Expression Region</b> | 2-489  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |