



Recombinant Escherichia coli P-protein (pheA)

Product Code	CSB-BP364342ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A9J8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTSENPLLAL REKISALDEK LLALLAERRE LAVEVGKAKL LSHRPVRDID RERDLLERLI TLGKAHHLDA HYITRLFQLI IEDSVLTQQA LLQQHLNKIN PHSARIAFLG PKGSYSHLAA RQYAARHFEQ FIESGCAKFA DIFNQVETGQ ADYAVVPIEN TSSGAINDVY DLLQHTSLSI VGEMTLTIDH CLLVSGTTDL STINTVYSHP QPFQQCSKFL NRYPHWKIEY TESTSAAMEK VAQAKSPHVA ALGSEAGGTL YGLQVLERIE ANQRQNFTRF VVLARKAINV SDQVPAKTTL LMATGQQAGA LVEALLVLRN HNLIMTRLES RPIHGNPWEE MFYLDIQANL ESAEMQKALK ELGEITRSMK VLGCPYSENV VPVDPT
Source	Baculovirus
Target Names	pheA
Protein Names	Recommended name: P-protein Including the following 2 domains: Chorismate mutase Short name= CM EC= 5.4.99.5 Prephenate dehydratase Short name= PDT EC= 4.2.1.51
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.