



# Recombinant Escherichia coli P-protein (pheA)

<b>Product Code</b>	CSB-MP364342ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A9J8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSENPLLAL REKISALDEK LLALLAERRE LAVEVGKAKL LSHRPVRDID RERDLLERLI TLGKAHHLDA HYITRLFQLI IEDSVLTQQA LLQQHLNKIN PHSARIAFLG PKGSYSHLAA RQYAARHFEQ FIESGCAKFA DIFNQVETGQ ADYAVVPIEN TSSGAINDVY DLLQHTSLSI VGEMTLTIDH CLLVSGTTDL STINTVYSHP QPFQQCSKFL NRYPHWKIEY TESTSAAMEK VAQAKSPHVA ALGSEAGGTL YGLQVLERIE ANQRQNFTRF VVLARKAINV SDQVPAKTTL LMATGQQAGA LVEALLVLRN HNLIMTRLES RPIHGPNWEE MFYLDIQANL ESAEMQKALK ELGEITRSMK VLGCYPSENV VPVDPT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pheA
<b>Protein Names</b>	Recommended name: P-protein Including the following 2 domains: Chorismate mutase Short name= CM EC= 5.4.99.5 Prephenate dehydratase Short name= PDT EC= 4.2.1.51
<b>Expression Region</b>	1-386
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.