



# Recombinant P-protein (pheA)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP364343SZB-B                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                             |
| <b>Uniprot No.</b>       | P0A9K0                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen Species</b> | Shigella flexneri                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Sequence</b>          | MTSENPLLAL REKISALDEK LLALLAERRE LAVEVGKAKL LSHRPVRDID<br>RERDLLERLI TLGKAHHLDA HYITRLFQLI IEDSVLTQQA LLQQHLNKKIN<br>PHSARIAFLG PKGSYSHLAA RQYAARHFEQ FIESGCAKFA DIFNQVETGQ<br>ADYAVVPIEN TSSGAINDVY DLLQHTSLSI VGEMTLTIDH CLLVSGTTDL<br>STINTVYSHP QPFQQCSKFL NRYPHWKIEY TESTSAAMEK VAQAKSPHVA<br>ALGSEAGGTL YGLQVLERIE ANQRQNFTRF VVLARKAINV SDQVPAKTTL<br>LMATGQQAGA LVEALLVLRN HNLIMTRLES RPIHGPNWEE MFYLDIQANL<br>ESAEMQKALK ELGEITRSMK VLGCPYSENV VPVDPT |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Names</b>      | pheA                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Names</b>     | Recommended name: P-protein Including the following 2 domains: Chorismate mutase Short name= CM EC= 5.4.99.5 Prephenate dehydratase Short name= PDT EC= 4.2.1.51                                                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-386                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                            |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                    |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                             |