



Recombinant Human Ubiquitin carboxyl-terminal hydrolase 2 (USP2)

Product Code	CSB-MP025710HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75604
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSQLSSTLKR YTESARYTDA HYAKSGYGAY TPSSYGANLA ASLLEKEKLG FKPVPTSSFL TRPRTYGPSS LLDYDRGRPL LRPDITGGGK RAESQTRGTE RPLGSGLSGG SGFPYGVTTN CLSYLPINAY DQGVTLTQKL DSQSDLARDF SSLRTSDSYR IDPRNLGRSP MLARTRKELC TLQGLYQTAS CPEYLVLDYLE NYGRKGSASQ VPSQAPPSRV PEIISPTYRP IGRYTLWETG KGQAPGPSRS SSPGRDGMNS KSAQGLAGLR NLGNTCFMNS ILQCLSNTRE LRDYCLQRLY MRDLHHGSNA HTALVEEFAK LIQTIWTSSP NDVVPSEFK TQIQRYAPRF VGYNQQDAQE FLRFLLDGLH NEVNRVTLRP KSNPENLDHL PDDEKGRQMW RKYLEREDSR IGDLVFGQLK SSLTCTDCGY CSTVDFPFWD LSLPIAKRGY PEVTLMDCMR LFTKEDVLDG DEKPTCCRCR GRKRCIKKFS IQRFPKILVL HLKRFSESRI RTSKLTTFVN FPLRDLDLRE FASENTNHAV YNLYAVSNHS GTTMGGHYTA YCRSPGTGEW HTFNDSSVTP MSSSQVRTSD AYLLFYELAS PPSRM
Source	Mammalian cell
Target Names	USP2
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase 2 EC= 3.4.19.12 Alternative name(s): 41 kDa ubiquitin-specific protease Deubiquitinating enzyme 2 Ubiquitin thioesterase 2 Ubiquitin-specific-processing protease 2
Expression Region	1-605
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.