



Recombinant Human GTPase Era, mitochondrial (ERAL1)

Product Code	CSB-EP007759HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75616
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	VSCVAGS AFSGPRLASA SRSNGQGSAL DHFLGFSQPD SSVTPCVPAV SMNRDEQDVL LVHHPDMPEN SRVLRVLLG APNAGKSTLS NQLLGRKVFP VSRKVHTTRC QALGVITEKE TQVILLDTPG IISPGKQKRH HLELSLLEDP WKSMEADLV VVLVDVSDKW TRNQLSPQLL RCLTKYSQIP SVLVMNKVDC LKQKSVLLEL TAALTEGVVN GKKLKMRAF HSHPGTHCPS PAVKDPNTQS VGNPQRIGWP HFKEIFMLSA LSQEDVKTLK QYLLTQAQPG PWEYHSAVLT SQTPEEICAN IIREKLEHL PQEVPYNVQQ KTAVWEEGPG GELVIQQKLL VPKESYVKLL IGPKGHVISQ IAQEAGHDLM DIFLCDVDIR LSVKLLK
Source	E.coli
Target Names	ERAL1
Protein Names	Recommended name: GTPase Era, mitochondrial Short name= H-ERA Short name= hERA Alternative name(s): Conserved ERA-like GTPase Short name= CEGA ERA-W ERA-like protein 1
Expression Region	44-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.