




Recombinant Human Cancer/testis antigen 2 (CTAG2)

Product Code	CSB-EP006127HU
Abbreviation	Recombinant Human CTAG2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75638
Alias	Autoimmunogenic cancer/testis antigen NY-ESO-2;Cancer/testis antigen 6.2 ;CT6.2L antigen family member 1 ;LAGE-1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	<pre>MQAEGQGTGGSTGDADGPGGPGIPDGPGGNAGGPGEAGATGGRGPRGAG AARASGPRGGAPRGPHGGAASAQDGRCPGARRPDSRLLQLHITMPFSSPM EAELVRRILSRDAAPLPRPGAVLKDFTVSGNLLFMSVRDQDREGAGRMRVVG WGLGSASPEGQKARDLRTPKHKVSEQRPGTPGPPPPEGAQGDGCRGVAFN VMFSAPHI</pre>
Research Area	Cancer
Source	E.coli
Target Names	CTAG2
Expression Region	1-210aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	37.1kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.