



Recombinant Human Ret finger protein-like 2 (RFPL2)

Product Code	CSB-MP019599HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75678
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEVAELGFPE TAVSQSRICL CAVLCGHWDF ADMMVIRSLS LIRLEGVEGR DPVGGGNLTN KRPSCAPSPQ DLSAQWKQLE DRGASSRRVD MAALFQEASS CPVCSDYLEK PMSLECGCAV CLKCINSLQK EPHGEDLLCC CSSMVSRRKNK IRRNRQLERL ASHIKELEPK LKKILQMNPR MRKFQVDMTL DANTANNFLL ISDDLRSVRS GRIRQNRQDL AERFDVSVCI LGSPRFTCGR HCWEVDVGTS TEWDLGVCRE SVHRKGRIQL TTELGFWTVS LRDGGRLSAT TVPLTFLFVD RKLQRVGIFL DMGMQNVSFF DAESGSHVYT FRSVSAEEPL RPFLAPSVPP NGDQGVLSIC PLMNSGTTDA PVRPGEAK
Source	Mammalian cell
Target Names	RFPL2
Protein Names	Recommended name: Ret finger protein-like 2 Alternative name(s): RING finger protein 79
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.