



Recombinant Escherichia coli Ribosomal large subunit pseudouridine synthase C (rluC)

Product Code	CSB-BP364430ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AA39
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKTETPSVKI VAITADEAGQ RIDNFLRTQL KGVPKSMIYR ILRKGEVRVN KKRIKPEYKL EAGDEVRIIP VRVAEREEEA VSPHLQKVAA LADVILYEDD HILVLNKPSG TAVHGGSGLS FGVIEGLRAL RPEARFLELV HRLDRDTSGV LLVAKKRSAL RSLHEQLREK GMQKDYALV RGQWQSHVKS VQAPLLKNIL QSGERIVRVS QEGKPSETRF KVEERYAFAT LVRCSPTGR THQIRVHTQY AGHPIAFDDR YGDREFDRQL TEAGTGLNRL FLHAAALKFT HPGTGEVMRI EAPMDEGLKR CLQKLRNAR
Source	Baculovirus
Target Names	rluC
Protein Names	Recommended name: Ribosomal large subunit pseudouridine synthase C EC=5.4.99.24 Alternative name(s): 23S rRNA pseudouridine(955/2504/2580) synthase rRNA pseudouridylate synthase C rRNA-uridine isomerase C
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.