



Recombinant 3-oxoacyl-[acyl-carrier-protein] synthase 2 (fabF)

Product Code	CSB-YP364501EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AAI6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	SKRRVVVTG LGMLSPVGNT VESTWKALLA GQSGISLIDH FDTSAYATKF AGLVKDFNCE DIISRKEQRK MDAFIQYGIV AGVQAMQDSG LEITEENATR IGAAIGSGIG GLGLIEENHT SLMNGGPRKI SPFFVPSTIV NMVAGHLTIM YGLRGPSISI ATRACTSGVHN IGHAARIAY GDADVMVAGG AEKASTPLGV GGFGAARALS TRNDNPQAAS RPWDKERDGF VLGDGAGMLV LEEYEHAKKR GAKIYAEVVG FGMSSDAYHM TSPPEENGAGA ALAMANALRD AGIEASQIGY VNAHGTSTPA GDKAEAQAVK TIFGEAASRV LVSSTKSMTG HLLGAAGAVE SIYSILALRD QAVPPTINLD NPDEGCDFD VPHEARQVSG MEYTLCSFG FGGTNGSLIF KKI
Source	Yeast
Target Names	fabF
Protein Names	Recommended name: 3-oxoacyl-[acyl-carrier-protein] synthase 2 EC= 2.3.1.179 Alternative name(s): 3-oxoacyl-[acyl-carrier-protein] synthase II Beta-ketoacyl-ACP synthase II Short name= KAS II
Expression Region	2-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.