



Recombinant Human Ras association domain-containing protein 9 (RASSF9)

Product Code	CSB-YP019381HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75901
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPFGRNLLK TRHKNRSPTK DMDSEEKEIV VWVCQEEKLV CGLTKRRTSA DVIQALLEEH EATFGEKRFL LGKPSDYCII EKWRGSERVL PPLTRILKLW KAWGDEQPNM QFVLVKADAF LPVPLWRTAE AKLVQNTEKL WELSPANMYK TLPPDKQKRI VRKTRKRLAK IKQDTVSHDR DNMETLVHLI ISQDHTIHQQ VKRMKELDLE IEKCEAKFHL DRVENDGENY VQDAYLMPSF SEVEQNLDLQ YEENQTLEDL SESDGIEQLE ERLKYRILI DKLSAEIEKE VKSVCIDINE DAEGEAASEL ESSNLESVKC DLEKSMKAGL KIHSHLSGIQ KEIKYSDSL QMKAKEYELL AKEFNLSLHIS NKDGCQLKEN RAKESEVPSS NGEIPPFTQR VFSNYTNDTD SDTGISSNHS QDSETTVGDV VLLST
Source	Yeast
Target Names	RASSF9
Protein Names	Recommended name: Ras association domain-containing protein 9 Alternative name(s): PAM COOH-terminal interactor protein 1 Short name= P-CIP1 Peptidylglycine alpha-amidating monooxygenase COOH-terminal interactor
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein localizes to perinuclear endosomes. This protein associates with peptidylglycine alpha-amidating monooxygenase, and may be involved with the trafficking of this enzyme through secretory or endosomal pathways.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.