



Recombinant Human Pre-mRNA-splicing factor SPF27 (BCAS2)

Product Code	CSB-BP002597HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75934
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AGTGLVAGE VVVDALPYFD QGYEAPGVRE AAAALVEEET RRYRPTKNYL SYLTAPDYS A FETDIMRNEF ERLAARQPIE LLSMKRYELP APSSGQKNDI TAWQECVNNS MAQLEHQAVR IENLELMSQH GCNAWKVYNE NLVHMIEHAQ KELQKLRKHI QDLNWQRKNM QLTAGSKLRE MESNWVSLVS KNYEIERTIV QLENEIYQIK QQHGEANKEN IRQDF
Source	Baculovirus
Target Names	BCAS2
Protein Names	Recommended name: Pre-mRNA-splicing factor SPF27 Alternative name(s): Breast carcinoma-amplified sequence 2 DNA amplified in mammary carcinoma 1 protein Spliceosome-associated protein SPF 27
Expression Region	2-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.