



Recombinant Uncharacterized protein ybhH (ybhH)

Product Code	CSB-MP364563EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AAV9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKKIPCVMMR GGTSRGAFLL AEHLPEDQTQ RDKILMAIMG SGNDLEIDGI GGGNPLTSKV AIISRSSDPR ADVDYLFQV IVHEQRVDTT PNCGNMLSGV GAFIENGLI AATSPVTRVR IRNVNTGTFI EADVQTPNGV VEYEGSARID GVPGTAAPVA LTFLNAAGTK TGKVFPTDNQ IDYFDDVPVT CIDMAMPVVI IPAAYLGKTG YELPAELDAD KALLARIESI RLQAGKAMGL GDVSNMVIPK PVLISPAQKG GAINVRYFMP HSCHRALAIT GAIAISSSCA LEGTVTRQIV PSVGYGNINI EHPSGALDVH LSNEGQDATT LRASVIRTTR KIFSGEVYLP
Source	Mammalian cell
Target Names	ybhH
Protein Names	Recommended name: Uncharacterized protein ybhH
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.