



Recombinant Human SAM and SH3 domain-containing protein 3 (SASH3)

Product Code	CSB-EP020715HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75995
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLRRKPSNAS EKEPTQKKKL SLQRSSSFKD FAKSKPSSPV VSEKEFNLDD NIPEDDSGVP TPEDAGKSGK KLGKKWRAVI SRTMNRKMGK MMVKALSEEM ADTLEEGSAS PTSPDYSLDS PGPEKMALAF SEQEEHELPV LSRQASTGSE LCSPSPGSGS FGEEPPAPQY TGPFCGRARV HTDFTPSPYD HDCLKLQKGD VIQIEKPPV GTWLGLLNGK VGSFKFIYVD VLPEEAVGHA RPSRRQSKGK RPKPKTLHEL LERIGLEEHT STLLNGYQT LEDFKELRET HLNELNIMDP QHRAKLLTAA ELLLDYDTGS EEAEEGAESS QEPVAHTVSE PKVDIPRDSG CFEGSESGRD DAELAGTEEQ LQGLSLAGAP
Source	E.coli
Target Names	SASH3
Protein Names	Recommended name: SAM and SH3 domain-containing protein 3 Alternative name(s): SH3 protein expressed in lymphocytes homolog
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.