



Recombinant Human Glutaredoxin-3 (GLRX3)

Product Code	CSB-YP009523HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O76003
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAGAAEAAVAAVEEVGSAGQFEELLRLKAKSLLVHFWAPWAPQCAQMNE VMAELAKELPQVSFVKLEAEGVPEVSEKYEISSVPTFLFFKNSQKIDRLDGAHA PELTKKVQRHASSGSFLPSANEHLKEDLNLRLKKLTHAAPCMLFMKGTQPQEP CGFSKQMVEILHKHNIQFSSFDIFSDEEVRQGLKAYSSWPTYYPQLYVSGELIGG LDIIKELEASEELDTICPKAPKLEERLKVLTNKASVMLFMKGNKQEAACGFSKQI LEILNSTGVEYETFDILEDEEVRQGLKAYSNWPTYYPQLYVKGELVGGGLDIVKEL KENGELLPILRGEN
Source	Yeast
Target Names	GLRX3
Protein Names	Recommended name: Glutaredoxin-3 Alternative name(s): PKC-interacting cousin of thioredoxin Short name= PICOT PKC-theta-interacting protein Short name= PKCq-interacting protein Thioredoxin-like protein 2
Expression Region	1-335aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.