



Recombinant Penicillin-binding protein activator LpoB (lpoB)

Product Code	CSB-MP364603EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AB39
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	C VGQREPAPVE EVKPAPEQPA EPQQPVPTVP SVPTIPQQPG PIEHEDQTAP PAPHIRHYDW NGAMQPMVSK MLGADGVTAG SVLLVDSVNN RTNGSLNAAE ATETLRNALA NNGKFTLVSA QQLSMAKQQL GLSPQDSLGT RSKAIGIARN VGAHYVLYSS ASGNVNAPTL QMQLMLVQTG EIIWSGKGAV SQQ
Source	Mammalian cell
Target Names	lpoB
Protein Names	Recommended name: Penicillin-binding protein activator LpoB Short name= PBP activator LpoB
Expression Region	20-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.