



Recombinant Uncharacterized protein yciM (yciM)

Product Code	CSB-EP364613EGX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AB59
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	WYMG RRSQQNKQD EANRLSRDYV AGVNFLLSNQ QDKAVDLFLD MLKEDTGTVE AHLTLGNLFR SRGEVDRAIR IHQTLMESAS LTYEQRLLAI QQLGRDYMAA GLYDRAEDMF NQLTDETDNR IGALQQLLQI YQATSEWQKA IDVAERLVKL GKDKQRVEIA HFYCELALQH MASDDLDRAM TLLKKGAAAD KNSARVSIMM GRVFMAGKEY AKAVESLQRV ISQDRELVSE TLEMLQTCYQ QLGKTAEWAE FLQRAVEENT GADAEMLLAD IIEARDGSEA AQVYITRQLQ RHPTMRVFK LMDYHLNEAE EGRAKESLMV LRDMVGEKVR SKPRYRCQKC GFTAYTLYWH CPSCRAWSTI KPIRGLDGL
Source	E.coli
Target Names	lapB
Protein Names	Recommended name: Uncharacterized protein yciM
Expression Region	17-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.