



# Recombinant Aspartate ammonia-lyase (aspA)

<b>Product Code</b>	CSB-EP364745SZB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0AC40
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Shigella flexneri
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSNNIRIEED LLGTREVPAD AYYGVHTLRA IENFYISNNK ISDIPEFVRG MVMVKKAAAM ANKELQTIPK SVANAIIAAC DEVLNNGKCM DQFPVDVYQG GAGTSVNMNT NEVLANIGLE LMGHQKGEYQ YLNPNDHVNK CQSTNDAYPT GFRIAVYSSL IKLVDAINQL REGFERKAVE FQDILKMGRT QLQDAVPMTL GQEFRAFSIL LKEEVKNIQR TAELLLEVNL GATAIGTGLN TPKEYSPLAV KKLAEVTGFP CVPAEDLIEA TSDCGAYVMV HGALKRLAVK MSKICNDLRL LSSGPRAGLN EINLPELQAG SSIMPAKVNP VVPEVVNQVC FKVIGNDTTV TMAAEAGQLQ LNVMEPVIGQ AMFESVHILT NACYNLLEKC INGITANKEV CEGYVYNSIG IVTYLNPFIG HHNGDIVGKI CAETGKSVRE VVLERGLLTE AELDDIFSVQ NLMHPAYKAK RYTDESEQ
<b>Source</b>	E.coli
<b>Target Names</b>	aspA
<b>Protein Names</b>	Recommended name: Aspartate ammonia-lyase Short name= Aspartase EC= 4.3.1.1
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.