



# Recombinant Aspartate ammonia-lyase (aspA)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP364745SZB  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | P0AC40   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Shigella flexneri  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MSNNIRIEED LLGTREVPAD AYYGVHTLRA IENFYISNNK ISDIPEFVRG<br>MVMVKKAAAM ANKELQTIPK SVANAIIAAC DEVLNNGKCM DQFPVDVYQG<br>GAGTSVNMNT NEVLANIGLE LMGHQKGEYQ YLNPNDHVNK CQSTNDAYPT<br>GFRIAVYSSL IKLVDAINQL REGFERKAVE FQDILKMGRT QLQDAVPMTL<br>GQEFRAFSIL LKEEVKNIQR TAELLLEVNL GATAIGTGLN TPKEYSPLAV<br>KKLAEVTGFP CVPAEDLIEA TSDCGAYVMV HGALKRLAVK MSKICNDLRL<br>LSSGPRAGLN EINLPELQAG SSIMPAKVNP VVPEVVNQVC FKVIGNDTTV<br>TMAAEAGQLQ LNVMEPVIGQ AMFESVHILT NACYNLLEKC INGITANKEV<br>CEGYVYNSIG IVTYLNPFIG HHNGDIVGKI CAETGKSVRE VVLERGLLTE<br>AELDDIFSVQ NLMHPAYKAK RYTDESEQ |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | aspA   |
| <b>Protein Names</b>     | Recommended name: Aspartate ammonia-lyase Short name= Aspartase EC= 4.3.1.1  |
| <b>Expression Region</b> | 1-478  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |