



Recombinant Escherichia coli Uncharacterized GST-like protein yibF (yibF)

Product Code	CSB-YP364767ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ACA1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKLVGSYTSP FVRKLSILLL EKGITFEFIN ELPYNADNGV AQFNPLGKVP VLVTEEGECW FDSPIIAEYI ELMNVAPAML PRDPLESLRV RKIEALADGI MDAGLVSVRE QARPAAQQSE DELLRQREKI NRSLDVLEGY LVDGTLKTD VNLATIAIAC AVGYLNFRRV APGWCVDRPH LVKLVENLFS RESFARTEPP KA
Source	Yeast
Target Names	yibF
Protein Names	Recommended name: Uncharacterized GST-like protein yibF
Expression Region	1-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.