



Recombinant Escherichia coli Right origin-binding protein (rob)

Product Code	CSB-BP364805ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ACI0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MDQAGIIRDL LIWLEGHLDQ PLSLDNVAAK AGYSKWHLQR MFKDVTGHAI GAYIRARRLS KSAVALRLTA RPILDIALQY RFDSQQTFTR AFKKQFAQTP ALYRRSPEWS AFGIRPPLRL GEFTMPEHKF VTLEDTPLIG VTQSYSCSLE QISDFRHEMR YQFWHDFLGN APTIPPVLYG LNETRPSQDK DDEQEVFYTT ALAQQADGY VLTGHPVMLQ GGEYVMFTYE GLGTGVQEFI LTVYGTCTMPM LNLTRRKGQD IERYYP AEDA KAGDRPINLR CELLPIRR
Source	Baculovirus
Target Names	rob
Protein Names	Recommended name: Right origin-binding protein
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.