



Recombinant Deoxyribose operon repressor (deoR)

Product Code	CSB-MP364818EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ACK6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	METRREERIG QLLQELKRSD KLHLKDAAAL LGVSEMTIRR DLNNHSAPVV LLGGYIVLEP RSASHYLLSD QKSRLVEEKR RAAKLAATLV EPDQTLFFDC GTTTPWIEA IDNEIPFTAV CYSLNTFLAL KEKPHCRAFL CGGEFHASNA IFKPIDFQQT LNNFCPDIAF YSAAGVHVS KATCFNLEEL PVKHWAMMSMA QKHVLVDHS KFGKVRPARG GDLKRFDIVV SDCCPEDEYV KYAQTQRIKL MY
Source	Mammalian cell
Target Names	deoR
Protein Names	Recommended name: Deoxyribose operon repressor
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.