



Recombinant HTH-type transcriptional repressor CytR (cytR)

Product Code	CSB-EP364834EGX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ACN8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKAKKQETAA TMKDVALKAK VSTATVSRAL MNPDKVSQAT RNRVEKAARE VGYPQPMGR NVKRNESRTI LVIVPDICDP FFSEIRGIE VTAANHGYSLV LIGDCAHQNQ QEKTFIDLII TKQIDGMLLL GSRLPFDASI EEQRNLPPMV MANEFAPELE LPTVHIDNLT AAFDAVNYLY EQGHKRIGCI AGPEEMPLCH YRLQGYVQAL RRCGIMVDPQ YIARGDFTFE AGSKAMQQLL DLPQPPTAVF CHSDVMALGA LSQAKRQGLK VPEDLSIIGF DNIDLTQFCD PPLTTIAQPR YEIGREAMLL LLDQMGGQHV GSGSRLMDCE LIIRGSTRAL P
Source	E.coli
Target Names	cytR
Protein Names	Recommended name: HTH-type transcriptional repressor CytR
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.