



# Recombinant Escherichia coli HTH-type transcriptional repressor PurR (purR)

<b>Product Code</b>	CSB-MP364839ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0ACP7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATIKDVAKR ANVSTTTVSH VINKTRFVAE ETRNAVWAAI KELHYSPSAV ARSLKVNHTK SIGLLATSSE AAYFAEIIIEA VEKNCFQKGY TLILGNAWNN LEKQRAYLSM MAQKRVDGLL VMCSEYPEPL LAMLEEYRHI PMVVMDWGEA KADFTDAVID NAFEGGYMAG RYLIERGHRE IGVIPGPLER NTGAGRLAGF MKAMEEAMIK VPESWIVQGD FEPESGYRAM QQILSQPHRP TAVFCGGDIM AMGALCAADE MGLRVPQDVS LIGYDNVRNA RYFTPALTTI HQPKDSLGET AFNMLLDRIV NKREEPQSIE VHPRLIERRS VADGPFDRDYR R
<b>Source</b>	Mammalian cell
<b>Target Names</b>	purR
<b>Protein Names</b>	Recommended name: HTH-type transcriptional repressor PurR Alternative name(s): Pur regulon repressor Purine nucleotide synthesis repressor
<b>Expression Region</b>	2-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.