



Recombinant Rat Epsin-1 (Epn1)

Product Code	CSB-MP007740RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88339
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MSTSSLRRQM KNIVHNYSEA EIKVREATSN DPWGPSSSLM SEIADLTYNV VAFSEIMSMI WKRLNDHGKN WRHVYKAMTL MEYLIKTGSE RVSQQCKENM YAVQTLKDFQ YVDRDGKDQG VNVREKAKQL VALLRDEDRL REERAHAKT KEKLAQTATA SSAAVGSGPP PEAEQAWPQS SGEEELQLQL ALAMSKEEAD QPPSCGPEDD VQLQLALSLS REEHDKEERI RRGDDLRLQM AIEESKRETG GKEESSLMDL ADVFTTPAPP QASDPWGGPA SVPTAVPVA AASDPWGAPA VPPAADPWGG AAPTPASGDP WRPAAPTGPS VDPWGGTPAP AAGEGPTSDP WGSADGGAPV SGPPSSDPWA PAPA FSDPWG GSPAKPSSNG TAVGGFDTEP DEFSDFDRLR TALPTSGSST GELELLAGEV PARSPGAFDM SGVGGSLAES VGSPPPAATP TPTPPTKTP ESFLGPNAAL VDLDSLVS RP GTPPPGAKAS NPFLPSGAPA TGPSVTNPFQ PAPPATLTLN QLRLSPVPPV PGAPPTYISP LGGGGLPPM MPPGPPAPNT NPFL
Source	Mammalian cell
Target Names	Epn1
Protein Names	Recommended name: Epsin-1 Alternative name(s): EPS-15-interacting protein 1
Expression Region	1-575
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.