



Recombinant AMP nucleosidase (amn)

Product Code	CSB-EP365070EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AE13
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	MNNKGSGLTTP AQALDKLDAL YEQSVVALRN AIGNYITSGE LPDENARKQG LFVYPSLTVT WDGSTTNPPK TRAFGRFTHA GSYTTTITRP TLFRSYLNEQ LTLTYQDYGA HISVQPSQHE IPYPYVIDGS ELTLDRSMSA GLTRYFPTTE LAQIGDETAD GIYHPTTEFSP LSHFDARRVD FSLARLRHYT GTPVEHFQPF VLFTNYTRYV DEFVRWGCSQ ILDPDSPYIA LSCAGGNWIT AETEAPEEAI SDLAWKKHQM PAWHLITADG QGITLVNIGV GPSNAKTICD HLAVLRPDVW LMIGHCGGLR ESQAIGDYVL AHAYLRDDHV LDAVLPPDIP IPSIAEVQRA LYDATKLVSG RPKGEEVKQRL RTGTVVTTDD RNWELRYSAS ALRFNLSRAV AIDMESATIA AQQYRFRVPY GTLLCVSDKP LHGEIKLPGQ ANRFYEGAIS EHLQIGIRAI DLLRAEGDRL HSRKLRTFNE PPFR
Source	E.coli
Target Names	amn
Protein Names	Recommended name: AMP nucleosidase EC= 3.2.2.4
Expression Region	1-484
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.