



Recombinant Mouse WAS/WASL-interacting protein family member 3 (Wipf3)

Product Code	CSB-BP026116MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C7L0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPVPPPPPPP PPPPPPPPPP LGAPPPPPPLG APPPPPPPGP PVSTDTPSLR KPKDLKGRSAL LADIQQGTRL RKVTQINDRS APQIEGSKGT SKEGGAAGSN ARGGNTPPAL GDLFAGGFPV LRPAGQRDVS GKGSGQGP GS RAPSRLPIK AISGPLPAPA SPRLGNASET HSSARPVPPR PSVPAPPPPT PPPPPPPPLP PASPIKAQLV SPPAPPTKVN PSVVPPPLPC APPLPPPPPT PPPLPPASAL SDKAVRPQLA PLHLPIPPP LPLLPPCGYP GLHSEPNSPA QEVREPPAPP PPPPPPPPP LPTYASCSSR TAVAPPPLPG ANNSGSETPP PLPPKSPSFQ TQKALPTPPG APGPQAILQK KRRGPGTSGG KLNPPPAPPA RSPTTELSSK SQQPGGQLRN GGQHAIDDFE SKFTFHSMED FPPPDEYKPC QKIYPSKVPR SRTPGSWLQA EAAGQSSDDI KSRNSQLSLK ALPVR
Source	Baculovirus
Target Names	Wipf3
Protein Names	Recommended name: WAS/WASL-interacting protein family member 3 Alternative name(s): Corticosteroids and regional expression protein 16 homolog
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.