



Recombinant Curli production assembly/transport component CsgG (csgG)

Product Code	CSB-EP365116EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEA3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	CLTAP PKEAARPTLM PRAQSYKDLT HLPAPTGKIF VSVYNIQDET GQFKPYPASN FSTAVPQSAT AMLVTALKDS RWFIPLERQG LQNLLNERKI IRAAQENGTV AINNRIPLQS LTAANIMVEG SIIGYESNVK SGGVGARYFG IGADTQYQLD QIAVNLRVVN VSTGEILSSV NTSKTILSYE VQAGVFRFID YQRLLEGEVG YTSNEPVMMLC LMSAIETGVI FLINDGIDRG LWDLQNKAEK QNDILVKYRH MSVPPES
Source	E.coli
Target Names	csgG
Protein Names	Recommended name: Curli production assembly/transport component CsgG
Expression Region	16-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.