



Recombinant Shigella flexneri Probable phosphoribulokinase (prkB)

Product Code	CSB-EP365215SZB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEX6
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	MSAKHPVIAV TGSSGAGTTT TSLAFRKIFA QLN LHAAEVE GDSFHRYTRP EMDMAIRKAR DAGRHISYFG PEANDFGLLE QTFIEYGQSG KGKSRKYLHT YDEAVPWNQV PGTFTPWQPL PEPTDVL FYE GLHGGVVTPQ HNVAQHVDLL VGVVPIVNLE WIQKLIRDTS ERGHSREAVM DSVVRSMEDY INYITPQFSR THLNFQRVPT VDTSNPFAAK GIPSLDES FV VIHFRNLEGI DFPWLLAMLQ GSFISHINTL VVPGGKMGLA MELIMLPLVQ RLMEGKKIE
Source	E.coli
Target Names	prkB
Protein Names	Recommended name: Probable phosphoribulokinase Short name= PRK Short name= PRKase EC= 2.7.1.19 Alternative name(s): Phosphopentokinase
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.