



Recombinant Escherichia coli O157:H7 Nucleoside triphosphate pyrophosphohydrolase (mazG)

Product Code	CSB-MP365219EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEY4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MNQIDRLLTI MQRLRDPENG CPWDKEQTFA TIAPYTTLEET YEVLDAIARE DFDDLRLGELG DLLFQVVFYA QMAQEEGRFD FNDICAAISD KLERRHPHVF ADSSAENSSE VLARWEQIKT EERAQKAQHS ALDDIPRSLP ALMRAQKIQK RCANVGFDWT TLGPVVDKVV EEIDEVMYEA RQAVVDQAKL EEEMGDLLFA TVNLARHLGT KAEIALQKAN EKFERRFREV ERIVAARGLE MTGVDLETME EVWQQVQRQE IDL
Source	Mammalian cell
Target Names	mazG
Protein Names	Recommended name: Nucleoside triphosphate pyrophosphohydrolase Short name= NTP-PPase EC= 3.6.1.8
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.