



# Recombinant Escherichia coli Protein mrp (mrp)

<b>Product Code</b>	CSB-EP365231ENV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0AF08
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNEQSQAKSP EALRAMVAGT LANFQHPTLK HNLTTLKALH HVAWMDDTLH VELVMPFVWH SAFEELKEQC SAELLRITGA KAIDWKLSHN IATLKRVKNQ PGINGVKNII AVSSGKGGVG KSSTAVNLAL ALAAEGAKVG ILDADIYGPS IPTMLGAENQ RPTSPDGTHM APIMSHGLAT NSIGYLVTTDD NAMVWRGPMA SKALMQMLQE TLWPDLDYLV LDMPPGTGDI QLTLAQNIPV TGAVVVTPPQ DIALIDAKKG IVMFEKVEVP VLGIVENMSV HICSNCGHHE PIFGTGGAEK LAEKYHTQLL GQMPLHISLR EDLDKGTPTV ISRPESEFTA IYRQLADRVA AQLYWQGEVI PGEISFRAV
<b>Source</b>	E.coli
<b>Target Names</b>	mrp
<b>Protein Names</b>	Recommended name: Protein mrp
<b>Expression Region</b>	1-369
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.