



Recombinant Escherichia coli Nitrate reductase molybdenum cofactor assembly chaperone NarJ (narJ)

Product Code	CSB-BP365240ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AF26
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIELVIVSRL LEYPDAALWQ HQQEMFEAIA ASKNLPKEDA HALGIFLRDL TTMDPLDAQA QYSELFDRGR ATSLLLFEHV HGESRDRGQA MVDLLAQYEQ HGLQLNSREL PDHLPLYLEY LAQLPQSEAV EGLKDIAPIL ALLSARLQQR ESRYAVLFDL LLKLANTAID SDKVAEKIAD EARDDTQPAL DAVWEEEQVK FFADKGC GDS AITAHQRRFA GAVAPQYLN I TTGGQH
Source	Baculovirus
Target Names	narJ
Protein Names	Recommended name: Nitrate reductase molybdenum cofactor assembly chaperone NarJ Alternative name(s): Redox enzyme maturation protein NarJ
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.