



Recombinant Serine/threonine-protein kinase-transforming protein raf (V-RAF)

Product Code	CSB-MP365443MKM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00532
Product Type	Recombinant Protein
Immunogen Species	Murine sarcoma virus 3611
Purity	>85% (SDS-PAGE)
Sequence	EKNKIRPRGQ RDSSYYWKME ASEVMLSTRI GSGSFGTVYK GKWHGDVAVK ILKVVDPTPE QLQAFRNEVA VLRKTRHVNI LLFMGYMTKD NLAIVTQWCE GSSLYKHLHV QETKFQMFQL IDIARQTAQG MDYLHAKNII HRDMKSNNIF LHEGLTVKIG DFGLATVKSR WSGSQQVEQP TGSVLWMAPE VIRMQDDNPF SFQSDVYSYG IVLYELMAGE LPYAHINNRD QIIFMVGRGY ASPDLSRLYK NCPKAIKRLV ADCVKKVKEE RPLFPQILSS IELLQHSLPK INRSAPEPSL HRAAHTEDIN ACTLTTSPRL PVF
Source	Mammalian cell
Target Names	V-RAF
Protein Names	Recommended name: Serine/threonine-protein kinase-transforming protein raf EC= 2.7.11.1
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.