



Recombinant *Saccharomyces cerevisiae* Phosphoglycerate kinase (PGK1)

Product Code	CSB-EP365449SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00560
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SLSSKLSVQ DLDLKDGRVF IRVDFNVPLD GKKITSNQRI VAALPTIKYV LEHHPRYVVL ASHLGRPNGE RNEKYS LAPV AKELQSL LGK DVTFLNDCVG PEVEAAVKAS APGSVILLEN LRYHIEEEGS RKVDGQKVKA SKEDVQKFRH ELSSLADVYI NDAFGTAHRA HSSMVGFDLP QRAAGFLLEK ELKYFGKALE NPTRPFLAIL GGAKVADKIQ LIDNLLDKVD SIIIGGGMAF TFKKVLENTE IGDSIFDKAG AEIVPKLMEK AKAKGVEVVL PVDFIIADAF SADANTKTVT DKEGIPAGWQ GLDNGPESRK LFAATVAKAK TIVWNGPPGV FEF EKFAAGT KALLDEVVKS SAAGNTVIIG GGDTATVAKK YGVTDKISHV STGGGASLEL LEGKELPGVA FLSEKK
Source	E.coli
Target Names	PGK1
Protein Names	Recommended name: Phosphoglycerate kinase EC= 2.7.2.3
Expression Region	2-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.