



Recombinant Bacillus licheniformis Subtilisin Carlsberg (subC)

Product Code	CSB-EP365470BQT-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00780
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Bacillus licheniformis
Purity	>85% (SDS-PAGE)
Sequence	AQTVP YGIPLIKADK VQAQGFKGAN VKVAVLDTGI QASHPDLNVV GGASFVAGEA YNTDGNGHGT HVAGTVAALD NTTGVLGVAP SVSLYAVKVL NSSGSGTYSG IVSGIEWATT NGMDVINMSL GGPSGSTAMK QAVDNAYARG VVVAAAGNS GSSGNTNTIG YPAKYDSVIA VGAVDSNSNR ASFSSVGAEL EVMAPGAGVY STYPTSTYAT LNGTSMASPH VAGAAALILS KHPNLSASQV RNRLSSTATY LGSSFYYGKG LINVEAAAQ
Source	E.coli
Target Names	apr
Protein Names	Recommended name: Subtilisin Carlsberg EC= 3.4.21.62
Expression Region	106-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.