



Recombinant Streptomyces tendae Alpha-amylase inhibitor HOE-467A

Product Code	CSB-MP365519SQB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01092
Product Type	Recombinant Protein
Immunogen Species	Streptomyces tendae
Purity	>85% (SDS-PAGE)
Sequence	DTTVSEPAPS CVTLYQSWRY SQADNGCAQT VTVKVVYEDD TEGLCYAVAP GQITTVGDGY IGSFHGARYL ARCL
Source	Mammalian cell
Protein Names	Recommended name: Alpha-amylase inhibitor HOE-467A Alternative name(s): Tendamistat
Expression Region	31-104
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.