



Recombinant Transforming protein Myb (V-MYB)

Product Code	CSB-MP365520ASQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01104
Product Type	Recombinant Protein
Immunogen Species	Avian myeloblastosis virus
Purity	≥85% (SDS-PAGE)
Sequence	NRTDVQCQHR WQKVLNPELN KGPWTKEEDQ RVIEHVQKYG PKRWSDIAKH LKGRIGKQCR ERWHNHLNPE VKKTSWTEEE DRIIYQAHKR LGNRWAEIAK LLPGRDNAV KNHWNSTMRR KVEQEGYPQE SSKAGPPSAT TGFQKSSHLM AFAHNPPAGP LPGAGQAPLG SDYPYHIAE PQNVPGQIPY PVALHINIIN VPQPAAAAIQ RHYTDEDPEK EKRIKELELL LMSTENELKG QQALPTQNHT ANYPGWHSTT VADNTRTSGD NAPVSCLGEH HHCTPSPVD HGCLPEESAS PARCMIVHQS NILDNVKNLL EFAETLQLID SFLNTSSNHE NLNLDNPALT STPVCGHKMS VTTPFHKDQT FTEYRKMHHG AV
Source	Mammalian cell
Target Names	V-MYB
Protein Names	Recommended name: Transforming protein Myb
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.