



Recombinant Viral myc transforming protein (MYC)

Product Code	CSB-EP365521AST-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01110
Product Type	Recombinant Protein
Immunogen Species	Avian myelocytomatosis virus MC29
Purity	>85% (SDS-PAGE)
Sequence	MPLSASLPSK NYDYDYDSVQ PYFYFEEEEEE NFYLAAQQRG SELQPPAPSE DIWKKFELLP MPPLSPSRRS SLAAASCFPS TADQLEMVTE LLGGDMVNQS FICDPDES FVKSIIIQDCM WSGFSAAKL EKVVSEKLAT YQASRQEGGP AAASRPGPPP SGPPPPAGP AASAGLYLHD LGAAAADCID PSVVFYPLS ERAPRAAPP ANPAALLGVD TPPTTSSDSE EEQEEDDEID VVTLAEANES ESSTESSTEA SEEHCKPHHS PLVLKRCHVN IHQHNYAAPP STKVEYPAAK RLKLDSEGRVL KQISNNRKCS SPRTL DSEEN DKRRTHNVLE RQRRNELKLR FFALRDQIPE VANNEKAPKV VILKKATEYV LSLQSDEHKL IAEKEQLRRR REQLKHNLEQ LRNSRA
Source	E.coli
Target Names	MYC
Protein Names	Recommended name: Viral myc transforming protein Short name= v-Myc
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.