



# Recombinant Human Extracellular matrix protein 2 (ECM2)

<b>Product Code</b>	CSB-BP007384HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94769
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>KNEEI<sup>PR</sup>KQR RKIYHRRLRK SSTSHK<sup>HRSN</sup> RQLGIQ<sup>QTTV</sup> FTPVARLP<sup>IV</sup>  NFDYSME<sup>EKF</sup> ESFSSF<sup>PGVE</sup> SSYNVLP<sup>GKK</sup> GHCLVK<sup>GITM</sup> YNKAVWS<sup>PEP</sup>  CTTCLC<sup>SDGR</sup> VLCDETM<sup>CHP</sup> QRCPQ<sup>TVIPE</sup> GECCP<sup>VCSAT</sup> VSYLL<sup>SGIA</sup>  LNDRNE<sup>FSGD</sup> SSEQRE<sup>PTNL</sup> LHKQL<sup>PPPQV</sup> GMDRIV<sup>RKEA</sup> LQSEED<sup>EEVK</sup>  EEDTE<sup>QKRET</sup> PESRN<sup>QGQLY</sup> SEGDS<sup>RGGDR</sup> KQRP<sup>GEERRL</sup>  AHQQ<sup>RQGRE</sup> EEEDE<sup>EEEGE</sup> EGEED<sup>EEDDEE</sup> DPVR<sup>GDMFRM</sup>  PSRS<sup>PLPAPP</sup> RGTLR<sup>LPSGC</sup> SLSY<sup>RTISCI</sup> NAMLT<sup>QIPPL</sup> TAPQ<sup>ITSLEL</sup>  TGNSI<sup>ASIPD</sup> EAFN<sup>GLPNLE</sup> RLDLS<sup>KNNIT</sup> SSGIG<sup>PKAFK</sup> LLK<sup>KLMRLNM</sup>  DGNN<sup>LIQIPS</sup> QLPST<sup>LEELK</sup> VNENN<sup>LQAID</sup> EESL<sup>SDLNQL</sup> VTLE<sup>ELEGNL</sup>  SEAN<sup>VNPLAF</sup> KPLK<sup>SLAYLR</sup> LGKN<sup>KFRIIP</sup> QGLP<sup>GSIEEL</sup> YLENN<sup>QIEEI</sup>  TEIC<sup>FNHTRK</sup> INVIV<sup>LRYNK</sup> IEEN<sup>RIAPLA</sup> WINQ<sup>ENLESI</sup> DLSY<sup>NKLYHV</sup>  PSYL<sup>PKSLLH</sup> LVLL<sup>GNQIER</sup> IPGY<sup>VFGHME</sup> PGLE<sup>YLYLSF</sup> NKLAD<sup>DGMDR</sup>  VSFY<sup>GAYHSL</sup> RELFL<sup>DHNDL</sup> KSIPP<sup>GIQEM</sup> KALH<sup>FRLN</sup> NKIR<sup>NILPEE</sup>  ICNA<sup>EEDDDS</sup> NLEHL<sup>HLENN</sup> YIKIR<sup>EIPSY</sup> TFSC<sup>IRSYSS</sup> IVLK<sup>PQNIK</sup></p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	ECM2
<b>Protein Names</b>	Recommended name: Extracellular matrix protein 2 Alternative name(s): Matrix glycoprotein SC1/ECM2
<b>Expression Region</b>	21-699
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Extracellular domain
<b>Target Details</b>	ECM2 encodes extracellular matrix protein 2, so named because it shares extensive similarity with known extracellular matrix proteins.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.