



Recombinant Human Interferon alpha-1/13 (IFNA1)

Product Code	CSB-BP365592HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01562
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	CDLPETH SLDNRRTLML LAQMSRISPS SCLMDRHDFG FPQEEFDGNQ FQKAPAI SVL HELIQQIFNL FTTKDSSAAW DEDLLDKFCT ELYQQQLNDLE ACVMQEERVG ETPLMNADSI LAVKKYFRR I TLYLTEKKYS PCAWEVVRAE IMRSLSLSTN LQERLRRKE
Source	Baculovirus
Target Names	IFNA1;
Protein Names	Recommended name: Interferon alpha-1/13 Short name= IFN-alpha-1/13 Alternative name(s): Interferon alpha-D Short name= LeIF D
Expression Region	24-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.