



Recombinant Human Actin-binding LIM protein 3 (ABLIM3)

Product Code	CSB-BP001109HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94929
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNTSIPYQQN PYNPRGSSNV IQCYRCGDTG KGEVVRVHNN HFHIRCFTCQ VCGCGLAQSG FFFKNQEYIC TQDYQQLYGT RCDSCRDFIT GEVISALGRT YHPKCFVCSL CRKPFPIGDK VTFSGKECVC QTCSQSMAS KPIKIRGPSH CAGCKEEIKH GQSLLALDKQ WHVSCFKCQT CSVILTGEYI SKDGVVPYCES DYHAQFGIKC ETCDRYISGR VLEAGGKHYH PTCARCVRCH QMFTEGEEMY LTGSEVWHPI CKQAARAEKK LKHRRTSETS ISPPGSSIGS PNRVICA KVD NEILNYKDLA ALPKVKSIYE VQRPDLISYE PHSRYMSDEM LERCGYGESL GTLSPYSQDI YENLDRQR ASSPGYIDSP TYSRQGMSP FSRSPHHYYR SGPESGRSSP YHSQLDVRSS TPTSYPAPKH FHIPAGDSNI YRKPPIYKRH GDLSTATKSK TSEDISQTSK YSPIYSPDPY YASESEYWTY HGSPKVPRAR RFSSGGEEDD FDRSMHKLQS GIGRLILKEE MKARSSSYAD PWTPPRSSTS SREALHTAGY EMSLNGSPRS HYLADSDPLI SKSASLPAYR RNLHRTPSA DLFHYDSMNA VNWGMREYKI YPYELLVTT RGRNRLPKDV DRTRLERHLS QEEFYQVFGM TISEFDRLAL WKRNELKKQA RLF
Source	Baculovirus
Target Names	ABLIM3
Protein Names	Recommended name: Actin-binding LIM protein 3 Short name= abLIM-3 Alternative name(s): Actin-binding LIM protein family member 3
Expression Region	1-683
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.