



# Recombinant Human Actin-binding LIM protein 3 (ABLIM3)

<b>Product Code</b>	CSB-EP001109HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94929
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNTSIPYQQN PYNPRGSSNV IQCYRCGDTC KGEVVRVHNN HFHIRCFTCC VCGCGLAQSG FFFKNQEYIC TQDYQQLYGT RCDSCRDFIT GEVISALGRT YHPKCFVCSL CRKPFPIGDK VTFSGKECVC QTCSQSMAS KPIKIRGSPH CAGCKEEIKH GQSLALDKQ WHVSCFKCQT CSVILTGEYI SKDGVVYCES DYHAQFGIKC ETCDRYISGR VLEAGGKHYH PTCARCVRCH QMFTEGEEMY LTGSEVWHPI CKQAARAEEK LKHRRTSETS ISPPGSSIGS PNRVICAQVD NEILNYKDLA ALPKVKSIYE VQRPDLISYE PHSRYMSDEM LERCGYGESL GTLSPYSQDI YENLDRQR ASSPGYIDSP TYSRQGMSP FSRSPHHYYR SGPESGRSSP YHSQLDVRSS TPTSQAPKH FHIPAGDSNI YRKPPIYKRH GDLSTATKSK TSEDISQTSK YSPIYSPDPY YASESEYWTY HGSPKVPRAR RFSSGGEEDD FDRSMHKLQS GIGRLILKEE MKARSSSYAD PWTPPRSSTS SREALHTAGY EMSLNGSPRS HYLADSDPLI SKSASLPAYR RNGLHRTPSA DLFHYSMNA VNWGMREYKI YPYELLVTT RGRNRLPKDV DRTRLERHLS QEEFYQVFGM TISEFDRLAL WKRNELKKQA RLF
<b>Source</b>	E.coli
<b>Target Names</b>	ABLIM3
<b>Protein Names</b>	Recommended name: Actin-binding LIM protein 3 Short name= abLIM-3 Alternative name(s): Actin-binding LIM protein family member 3
<b>Expression Region</b>	1-683
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.