



Recombinant Hevea brasiliensis Pro-hevein (HEV1)

| | |
|--------------------------|---|
| Product Code | CSB-EP365797HWI-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P02877 |
| Product Type | Recombinant Protein |
| Immunogen Species | Hevea brasiliensis (Para rubber tree) (Siphonia brasiliensis) |
| Purity | >85% (SDS-PAGE) |
| Sequence | EQC GRQAGGKLC P NNLCCSQWGW CGSTDEYCSP DHNCQSNCKD SGEGVGGGSA SNVLATYHLY NSQDHGWDLN AASAYCSTWD ANKPYSWRSK YGWTAF CGPV GAHGQSSCGK CLSVTNTGTG AKTTVRIVDQ CSNGGLDL DV NVFRQLD TDG KGYERGHITV NYQFVDCGDS FNPLFSVMKS SVIN |
| Source | E.coli |
| Target Names | HEV1 |
| Protein Names | Recommended name: Pro-hevein Alternative name(s): Major hevein Cleaved into the following 2 chains: 1. Hevein Alternative name(s): Allergen= Hev b 6 Win-like protein |
| Expression Region | 18-204 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |