



Recombinant PI protein (pir)

Product Code	CSB-MP365836ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03067
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MRLKVMMDVN KGTKIRHRNE LNHTLAQLPL PAKRVMYMAL APIDSKEPLE RGRVFKIRAE DLAALAKITP SLAYRQLKEG GKLLGASKIS LRGDDIIALA KELNLPFTAK NSPEELDLNI IEWIAYSPDE GYLSLKFTRT IEPYISSLIG KKNKFTTQLL TASLRLSSQY SSSLYQLIRK HYSNFKKKNY FIISVDELKE ELIAYTFDKD GNIEYKYPDF PIFKRDVLNK AIAEIKKKTE ISFVGFTVHE KEGRKISKLK FEFVDEDEF SGDKDDEAFF MNLSEADAAF LKVFDVETVPP KKAKG
Source	Mammalian cell
Target Names	pir
Protein Names	Recommended name: PI protein Alternative name(s): Replication initiation protein
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.