



Recombinant Rous sarcoma virus Gag-Pro polyprotein (gag-pro)

Product Code	CSB-BP365907RHO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03322
Product Type	Recombinant Protein
Immunogen Species	Rous sarcoma virus (strain Prague C) (RSV-PrC)
Purity	>85% (SDS-PAGE)
Sequence	MEAVIKVISS ACKTYCGKTS PSKKEIGAML SLLQKEGLLM SPSPDLYSPGS WDPITAALSQ RAMILGKSGE LKTWGLVLGA LKAAREEQVT SEQAKFWLGL GGGRVSPPGP ECIEKPATER RIDKGEEVGE TTVQRDAKMA PEETATPKTV GTSCY
Source	Baculovirus
Target Names	gag
Protein Names	Recommended name: Gag-Pro polyprotein Cleaved into the following 8 chains: 1. p3 2. Matrix protein p19 3. p2A 4. p2B 5. p10 6. Capsid protein p27 7. Nucleocapsid protein p12 8. Protease p15 EC= 9. 3.4.23.-
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.