



Recombinant Squirrel monkey retrovirus Pol polyprotein (pol), partial

| | |
|--------------------------|---|
| Product Code | CSB-BP365919SKS |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P03364 |
| Product Type | Recombinant Protein |
| Immunogen Species | Squirrel monkey retrovirus (SMRV-H) (SMRV-HLB) |
| Purity | >85% (SDS-PAGE) |
| Source | Baculovirus |
| Target Names | pol |
| Protein Names | Recommended name: Pol polyprotein Cleaved into the following 2 chains: 1. Reverse transcriptase/ribonuclease H Short name= 2. RT EC= 3. 2.7.7.49 EC= 4. 2.7.7.7 EC= 5. 3.1.26.4 6. Integrase Short name= 7. IN |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |