



Recombinant Enterobacteria phage alpha3 Minor spike protein (H)

Product Code	CSB-EP366008ECD-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03650
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage alpha3 (Bacteriophage alpha-3)
Purity	>85% (SDS-PAGE)
Sequence	MLGAVVGGIA SALASGAASK LFGGSQRGVS VMQQGAESAG LTNGGGAIMS DNDQGIQSAI QGSNVPPAGQ LPESQNSGVM ADAGNMLRNA GKSLLDGTIQ AGSDKVKQAL IDKLIGGNDK KDRGKATRDY LAAAFPELNP WERAGAGAST AGLESSAQNQ QKEMLMQLD NQKDIAKMQM DNNLQIAGIQ SATSRQNTKD SVYAQNEMLQ YNQRESQARV ASILANTDLT TKQATHEIMR MALTRAQETG QHLTNSQIMA LEKKVYAEIG KIHQDTQNSR YGSSQVTAAS KDVNTMITDA ASGASDWVSQ QWNSFFKDGK SNGISLNTRK
Source	E.coli
Target Names	H
Protein Names	Recommended name: Minor spike protein Alternative name(s): H protein Pilot protein
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.