



Recombinant Human Cellular communication network factor 6 (CCN6)

Product Code	CSB-BP026121HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95389
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TGPLDTT PEGRPGEVSD APQRKQFCHW PCKCPQQKPR CPPGVSLVRD GCGCKICAK QPGEICNEAD LCDPHKGLYC DYSVDRPRYE TGVCAYLVAV GCEFNQVHYH NGQVFQPNPL FSCLCVSGAI GCTPLFIPKL AGSHCSGAKG GKKSQSNCS LEPLLQQLST SYKTMPAYRN LPLIWKKKCL VQATKWTPCS RTCGMGISNR VTNENSNCCEM RKEKRLCYIQ PCDSNILKTI KIPKGKTCQP TFQLSKAEKF VFSGCSSTQS YKPTFCGICL DKRCCIPNKS KMITIQFDPCP NEGSFKWKML WITSCVCQRN CREPGDIFSE LKIL
Source	Baculovirus
Target Names	CCN6
Protein Names	Recommended name: WNT1-inducible-signaling pathway protein 3 Short name= WISP-3Alternative name(s): CCN family member 6
Expression Region	24-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.