



Recombinant Human Tumor necrosis factor receptor superfamily member 6B (TNFRSF6B)

Product Code	CSB-YP023982HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95407
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	V AETPTYPWDR AETGERLVCA QCPPGTFVQR PCRRDSPTTC GPCPPRHYTQ FWNYLERCRY CNVLCGEREE EARACHATHN RACRCRTGFF AHAGFCLEHA SCPPGAGVIA PGTPSQNTQC QPCPPGTFSA SSSSSEQCQP HRNCTALGLA LNVPGSSSHD TLCTSCTGFP LSTRVPGAEE CERAVIDFVA FQDISIKRLQ RLLQALEAPE GWGPTPRAGR AALQLKLRRR LTELLGAQDG ALLVRLQLAL RVARMPGLER SVRERFLPVH
Source	Yeast
Target Names	TNFRSF6B
Protein Names	Recommended name: Tumor necrosis factor receptor superfamily member 6B Alternative name(s): Decoy receptor 3 Short name= DcR3 Decoy receptor for Fas ligand M68
Expression Region	30-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene belongs to the tumor necrosis factor receptor superfamily. The encoded protein is postulated to play a regulatory role in suppressing FasL- and LIGHT-mediated cell death. It acts as a decoy receptor that competes with death receptors for ligand binding. Overexpression of this gene has been noted in gastrointestinal tract tumors, and it is located in a gene-rich cluster on chromosome 20, with other potentially tumor-related genes. Two transcript variants encoding the same isoform, but differing in the 5' UTR, have been observed for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.